

Dilation

CRE™ Balloon Dilatation Catheter

- CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters..... 57
- CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters..... 57
- CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters..... 58
- CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters..... 58
- CRE Fixed Wire Esophageal Balloon Dilatation Catheters 59

MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter

- MaxForce TTS Esophageal High Performance Balloon Dilatation Catheters..... 60

Rigiflex™ II Single-Use Achalasia Balloon

- Rigiflex II Single-Use Achalasia Balloon 60

Inflation System


- Alliance™ II Integrated Inflation / Lithotripsy Device..... 61
- Alliance II Single-Use Syringe / Gauge Assembly 61
- CRE SteriFlate Disposable Inflation Device..... 61

Dilation Accessories

- Inflation / Deflation Devices and Accessories 61

Dilation

CRE™ PRO Balloon Dilatation Catheters



CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters



| Order Number | GTIN | Balloon Inflated O.D. | | Balloon Length (cm) | Inflation Pressures | | Minimum Working Channel (mm) | Catheter Size (Fr) | Working Length (cm) | Packaging |
|--------------|----------------|-----------------------|------------|---------------------|---------------------|--------------|------------------------------|--------------------|---------------------|-----------|
| | | (mm) | (Fr) | | (atm) | (kPa) | | | | |
| M00558600 | 08714729797487 | 6-7-8 | 18-21-24 | 5.5 | 3-6-10 | 304-608-1013 | 2.8 / 3.2* | 7.5 | 180 | Each |
| M00558610 | 08714729797494 | 8-9-10 | 24-27-30 | 5.5 | 3-5.5-9 | 304-557-912 | 2.8 / 3.2* | 7.5 | 180 | Each |
| M00558620 | 08714729797500 | 10-11-12 | 30-33-36 | 5.5 | 3-5-8 | 304-507-811 | 2.8 / 3.2* | 7.5 | 180 | Each |
| M00558630 | 08714729797517 | 12-13.5-15 | 36-40.5-45 | 5.5 | 3-4.5-8 | 304-456-811 | 2.8 / 3.2* | 7.5 | 180 | Each |
| M00558640 | 08714729797524 | 15-16.5-18 | 45-49.5-54 | 5.5 | 3-4.5-7 | 304-456-709 | 2.8 / 3.7* | 7.5 | 180 | Each |
| M00558650** | 08714729797531 | 18-19-20 | 54-57-60 | 5.5 | 3-4.5-6 | 304-456-608 | 2.8 / 3.7* | 7.5 | 180 | Each |

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Syringe / Gauge Assembly Order #5060,

[\(Click here to view page 61\).](#)

Guidewire included and preloaded in balloon dilatation catheter.

*Duodenoscope channel size.

**Esophageal and Biliary application only.



CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters



| Order Number | GTIN | Balloon Inflated O.D. | | Balloon Length (cm) | Inflation Pressures | | Minimum Working Channel (mm) | Catheter Size (Fr) | Working Length (cm) | Packaging |
|--------------|----------------|-----------------------|------------|---------------------|---------------------|--------------|------------------------------|--------------------|---------------------|-----------|
| | | (mm) | (Fr) | | (atm) | (kPa) | | | | |
| M00558660 | 08714729797548 | 6-7-8 | 18-21-24 | 5.5 | 3-6-10 | 304-608-1013 | 2.8 / 3.2* | 7.5 | 240 | Each |
| M00558670 | 08714729797555 | 8-9-10 | 24-27-30 | 5.5 | 3-5.5-9 | 304-557-912 | 2.8 / 3.2* | 7.5 | 240 | Each |
| M00558680 | 08714729797562 | 10-11-12 | 30-33-36 | 5.5 | 3-5-8 | 304-507-811 | 2.8 / 3.2* | 7.5 | 240 | Each |
| M00558690 | 08714729797579 | 12-13.5-15 | 36-40.5-45 | 5.5 | 3-4.5-8 | 304-456-811 | 2.8 / 3.2* | 7.5 | 240 | Each |
| M00558700 | 08714729797586 | 15-16.5-18 | 45-49.5-54 | 5.5 | 3-4.5-7 | 304-456-709 | 2.8 / 3.7* | 7.5 | 240 | Each |
| M00558710** | 08714729797593 | 18-19-20 | 54-57-60 | 5.5 | 3-4.5-6 | 304-456-608 | 2.8 / 3.7* | 7.5 | 240 | Each |

Recommended Inflation Device: Alliance II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Syringe / Gauge Assembly Order #5060,

[\(Click here to view page 61\).](#)


Guidewire included and preloaded in balloon dilatation catheter.

*Duodenoscope channel size.

**Esophageal, Colonic and Biliary application only.

Dilation

CRE™ Balloon Dilatation Catheters



CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters



| Order Number | GTIN | Balloon Inflated O.D. | | Balloon Length (cm) | Inflation Pressures | | Minimum Working Channel (mm) | Catheter Size (Fr) | Working Length (cm) | Packaging |
|--------------|----------------|-----------------------|------------|---------------------|---------------------|--------------|------------------------------|--------------------|---------------------|-----------|
| | | (mm) | (Fr) | | (atm) | (kPa) | | | | |
| M00558390 | 08714729339328 | 6-7-8 | 18-21-24 | 5.5 | 3-6-10 | 304-608-1013 | 2.8 | 7.5 | 180 | Each |
| M00558400 | 08714729339335 | 8-9-10 | 24-27-30 | 5.5 | 3-5.5-9 | 304-557-912 | 2.8 | 7.5 | 180 | Each |
| M00558410 | 08714729339342 | 10-11-12 | 30-33-36 | 5.5 | 3-5-8 | 304-507-811 | 2.8 | 7.5 | 180 | Each |
| M00558420 | 08714729339359 | 12-13.5-15 | 36-40.5-45 | 5.5 | 3-4.5-8 | 304-456-811 | 2.8 | 7.5 | 180 | Each |
| M00558430 | 08714729339366 | 15-16.5-18 | 45-49.5-54 | 5.5 | 3-4.5-7 | 304-456-709 | 2.8 | 7.5 | 180 | Each |
| M00558440* | 08714729339373 | 18-19-20 | 54-57-60 | 5.5 | 3-4.5-6 | 304-456-608 | 2.8 | 7.5 | 180 | Each |

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, [\(Click here to view page 61\)](#).

Guidewire included and preloaded in balloon dilatation catheter.

*Esophageal application only.

CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters



| Order Number | GTIN | Balloon Inflated O.D. | | Balloon Length (cm) | Inflation Pressures | | Minimum Working Channel (mm) | Catheter Size (Fr) | Working Length (cm) | Packaging |
|--------------|----------------|-----------------------|------------|---------------------|---------------------|--------------|------------------------------|--------------------|---------------------|-----------|
| | | (mm) | (Fr) | | (atm) | (kPa) | | | | |
| M00558450 | 08714729339380 | 6-7-8 | 18-21-24 | 5.5 | 3-6-10 | 304-608-1013 | 2.8 | 7.5 | 240 | Each |
| M00558460 | 08714729339397 | 8-9-10 | 24-27-30 | 5.5 | 3-5.5-9 | 304-557-912 | 2.8 | 7.5 | 240 | Each |
| M00558470 | 08714729339403 | 10-11-12 | 30-33-36 | 5.5 | 3-5-8 | 304-507-811 | 2.8 | 7.5 | 240 | Each |
| M00558480 | 08714729339410 | 12-13.5-15 | 36-40.5-45 | 5.5 | 3-4.5-8 | 304-456-811 | 2.8 | 7.5 | 240 | Each |
| M00558490 | 08714729339427 | 15-16.5-18 | 45-49.5-54 | 5.5 | 3-4.5-7 | 304-456-709 | 2.8 | 7.5 | 240 | Each |
| M00558500* | 08714729339434 | 18-19-20 | 54-57-60 | 5.5 | 3-4.5-6 | 304-456-608 | 2.8 | 7.5 | 240 | Each |


Recommended Inflation Device: Alliance II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, [\(Click here to view page 61\)](#).

Guidewire included and preloaded in balloon dilatation catheter.

*Esophageal and Colonic application only.

Dilation

CRE™ Balloon Dilatation Catheters (Continued)



CRE Fixed Wire Esophageal Balloon Dilatation Catheters




| Order Number | GTIN | Balloon Inflated O.D. | | Balloon Length (cm) | Inflation Pressures | | Minimum Working Channel (mm) | Catheter Size (Fr) | Working Length (cm) | Packaging |
|--------------|----------------|-----------------------|------------|---------------------|---------------------|--------------|------------------------------|--------------------|---------------------|-----------|
| | | (mm) | (Fr) | | (atm) | (kPa) | | | | |
| M00558330 | 08714729201984 | 6-7-8 | 18-21-24 | 8 | 3-6-10 | 304-608-1013 | 2.8 | 6 | 180 | Each |
| M00558340 | 08714729201991 | 8-9-10 | 24-27-30 | 8 | 3-5.5-9 | 304-557-912 | 2.8 | 6 | 180 | Each |
| M00558350 | 08714729202004 | 10-11-12 | 30-33-36 | 8 | 3-5-8 | 304-507-811 | 2.8 | 6 | 180 | Each |
| M00558360 | 08714729195986 | 12-13.5-15 | 36-40.5-45 | 8 | 3-4.5-8 | 304-456-811 | 2.8 | 6 | 180 | Each |
| M00558370 | 08714729195993 | 15-16.5-18 | 45-49.5-54 | 8 | 3-4.5-7 | 304-456-709 | 2.8 | 6 | 180 | Each |
| M00558380 | 08714729202011 | 18-19-20 | 54-57-60 | 8 | 3-4.5-6 | 304-456-608 | 2.8 | 6 | 180 | Each |
| M00558361 | 08714729339298 | 12-13.5-15 | 36-40.5-45 | 8 | 3-4.5-8 | 304-456-811 | 2.8 | 6 | 180 | Box 5 |
| M00558371 | 08714729339304 | 15-16.5-18 | 45-49.5-54 | 8 | 3-4.5-7 | 304-456-709 | 2.8 | 6 | 180 | Box 5 |
| M00558381 | 08714729339311 | 18-19-20 | 54-57-60 | 8 | 3-4.5-6 | 304-456-608 | 2.8 | 6 | 180 | Box 5 |

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060,

[▶ \(Click here to view page 61\).](#)

Dilation

MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter



MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter




| Order Number | GTIN | Balloon Inflated O.D. | | Balloon Length (cm) | Inflation Pressures | | Minimum Working Channel (mm) | Catheter Size | | Working Length (cm) | Packaging |
|--------------|----------------|-----------------------|------|---------------------|---------------------|-------|------------------------------|---------------|------|---------------------|-----------|
| | | (mm) | (Fr) | | (atm) | (kPa) | | (mm) | (Fr) | | |
| M00589060 | 08714729099581 | 6 | 18 | 6 | 10 | 1013 | 2.8 | 2.0 | 6 | 180 | Each |
| M00589080 | 08714729099598 | 8 | 24 | 6 | 10 | 1013 | 2.8 | 2.0 | 6 | 180 | Each |
| M00589100 | 08714729099895 | 10 | 30 | 6 | 7 | 709 | 2.8 | 2.0 | 6 | 180 | Each |
| M00589120 | 08714729099901 | 12 | 36 | 6 | 7 | 709 | 2.8 | 2.0 | 6 | 180 | Each |
| M00589140 | 08714729080602 | 14 | 42 | 6 | 7 | 709 | 2.8 | 2.0 | 6 | 180 | Each |
| M00589150 | 08714729099918 | 15 | 45 | 6 | 5 | 507 | 2.8 | 2.0 | 6 | 180 | Each |
| M00589160 | 08714729080619 | 16 | 48 | 6 | 5 | 507 | 2.8 | 2.0 | 6 | 180 | Each |
| M00589180 | 08714729099925 | 18 | 54 | 6 | 5 | 507 | 2.8 | 2.0 | 6 | 180 | Each |

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Sterile Syringe Order #5060,

[▶ \(Click here to view page 61\).](#)

Rigiflex™ II Single-Use Achalasia Balloon Dilatation Catheter



Rigiflex II Single-Use Achalasia Balloon Dilators*



| Order Number | GTIN | Balloon Inflated O.D. | | Balloon Length (cm) | Inflation Pressures | | Working Length (cm) | Catheter Size | | Recommended Guidewire** (in) | Packaging |
|--------------|----------------|-----------------------|------|---------------------|---------------------|-------|---------------------|---------------|------|------------------------------|-----------|
| | | (mm) | (Fr) | | (psi) | (kPa) | | (mm) | (Fr) | | |
| M00554500 | 08714729719106 | 30 | 90 | 10 | 20 | 138 | 90 | 4.6 | 14 | .038 | Each |
| M00554510 | 08714729719113 | 35 | 105 | 10 | 20 | 138 | 90 | 4.6 | 14 | .038 | Each |
| M00554520 | 08714729719120 | 40 | 120 | 10 | 20 | 138 | 90 | 4.6 | 14 | .038 | Each |


*Product Shipped Sterile. Recommended use with Achalasia Pneumatic Hand Pump and Monitor Order #5320,

[▶ \(Click here to view page 61\).](#)

**Guidewire included with Rigiflex II Balloon Dilatation Catheter.

Dilation

Inflation System



Alliance™ II Single-Use Syringe / Gauge Assembly / Integrated Inflation / Lithotripsy Device

| Order Number | GTIN | Description | Packaging |
|--------------|----------------|---|-----------|
| M00550620 | 08714729283881 | Alliance II Integrated Inflation / Lithotripsy Device | Each |
| M00550601 | 08714729129332 |  Alliance II Single-Use Syringe / Gauge Assembly | Box 5 |

Dilation Accessories

Inflation / Deflation Devices and Accessories

| Order Number | GTIN | Description | Packaging |
|--------------|----------------|---|-----------|
| M00553200 | 08714729195450 | Pneumatic Inflator | Each |
| M00550630 | 00856877002506 | CRE™ SteriFlate Disposable Inflation Device | Each |

Notes and Future Product Updates
