

Radial Jaw 4® Pulmonary

Single-Use Biopsy Forceps

Boston
Scientific

made with
surgical stainless steel

Samples
like never before™



Endoscopy

Introducing

Radial Jaw[®] 4 Single-Use Pulmonary Biopsy Forceps

The new Pulmonary Forceps now available in our family of market-leading single-use biopsy forceps



NEW SURGICAL STAINLESS STEEL JAW WITH IMPROVED MICROMESH TEETH

Designed to Provide:

- ▶ Tissue specimens for excellent sample handling and preparation
- ▶ Clean, precise bite for accurate histological diagnosis

NEW STREAMLINED CATHETER

Designed to Provide:

- ▶ Enhanced passability through tortuous anatomy
- ▶ The right balance of columnar strength and flexibility for excellent pushability and control during scope passage

DUAL FENESTRATION HOLE

- ▶ Designed to provide consistent sample retention and removal (Standard Capacity only)

SINGLE-USE

- ▶ Eliminates the risk of transmitting patient-to-patient disease
- ▶ Provides first time sharpness



NEW DISTAL END TUBE

- ▶ Improved visibility
- ▶ Prevents inadvertent lodging of the cap in the scope working channel



Boston Scientific, setting new standards for
pulmonary biopsy.

ORDERING INFORMATION

Order Number	Product Description	Differentiator	Jaw OD	Working Length (cm)	Minimum Working Channel	Box Quantity
M00515181	Single Use Biopsy Forceps	Pulmonary Standard Capacity	1.8mm	100cm	2.0mm	5
M00515182	Single Use Biopsy Forceps	Pulmonary Standard Capacity	1.8mm	100cm	2.0mm	20
M00515191	Single Use Biopsy Forceps	Pulmonary Standard Capacity w/ Needle	1.8mm	100cm	2.0mm	5
M00515192	Single Use Biopsy Forceps	Pulmonary Standard Capacity w/ Needle	1.8mm	100cm	2.0mm	20
M00515201	Single Use Biopsy Forceps	Pulmonary Large Capacity	2.2mm	100cm	2.8mm	5
M00515202	Single Use Biopsy Forceps	Pulmonary Large Capacity	2.2mm	100cm	2.8mm	20

Radial Jaw is a registered trademark of Boston Scientific Corporation or its affiliates.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Boston Scientific

Defining tomorrow, today.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com/endoscopy

Ordering Information
1.800.225.3226

© 2012 Boston Scientific Corporation or its affiliates. All rights reserved.

ENDO-88602-AA June 2012