


Boston
Scientific

ExtractorTM Pro

XL and DL Retrieval Balloons




*Defining tomorrow, today
in Endoscopy.*

Extractor™ Pro

XL and DL Retrieval Balloons

The Extractor Pro XL and DL Retrieval Balloons have been redesigned to offer long-wire users additional device options and a stiffer catheter for improved “pushability.”



Enhanced Balloon Design

- New and improved balloon is designed to improve strength and reliability.
- Each of the three balloon sizes can be inflated to two distinct diameters: 9–12mm; 12–15mm; 15–18mm to better accommodate anatomical variations and help obviate the need to use an additional balloon and change catheters if an alternative balloon diameter is required.
- “Squared shoulder” balloon is designed to facilitate stone removal.

Improved Catheter Design

- Catheter is designed to be easily back loaded over a guidewire.
- XL catheter enables front or back loading of guidewire and offers injection-above or injection-below options to support physician preference and facilitate procedural needs.
- Designed for excellent injection flow rate.



Ordering Information

Extractor Pro XL Triple-lumen Retrieval Balloon

Order Number	Balloon O.D. (mm)	Injection	Catheter O.D. (F/mm)
M00547100	9–12	Above	7–6 / 2.33–2
M00547110	12–15	Above	7–6 / 2.33–2
M00547120	15–18	Above	7–6 / 2.33–2
M00547130	9–12	Below	7–6 / 2.33–2
M00547140	12–15	Below	7–6 / 2.33–2
M00547150	15–18	Below	7–6 / 2.33–2

Packaged one per box

Recommended guidewire: .035" (0.89mm) Dreamwire™ High-Performance Guidewire. System length is 200cm.

Extractor Pro DL Double-lumen Retrieval Balloon

Order Number	Balloon O.D. (mm)	Injection	Catheter O.D. (F/mm)
M00547200	9–12	N/A	7–6 / 2.33–2
M00547210	12–15	N/A	7–6 / 2.33–2
M00547220	15–18	N/A	7–6 / 2.33–2

Packaged one per box

European Distribution Centre – The Netherlands
T: +31 45 54 67 700 F: +31 45 54 67 800
European Headquarters – Paris
T: +33 1 57 66 80 00 F: +33 1 57 66 84 99
Argentina (Freephone)
T: 0800 555 2678 F: +5411 4896 8517
Australia/New Zealand
T: +61 2 8063 8100 F: +61 2 9330 1404
Austria
T: +43 1 60 810 F: +43 1 60 810 60
Belgium (Freephone)
T: 0800 94 494 F: 0800 93 343
Brazil
T: +55 11 5502 8500 F: +55 11 5103 2212
Canada
T: +1 888 359 9691 F: +1 888 575 73965
Chile
T: +562 445 4904 F: +562 445 4915
China – Beijing
T: +86 10 8525 1588 F: +86 10 8525 1566
China – Guangzhou
T: +86 20 8767 9791 F: +86 20 8767 9789
China – Shanghai
T: +86 21 6141 5959 F: +86 21 6141 5900

Colombia
T: +57 1 629 5045 F: +57 1 612 4761
Czech Republic
T: +420 296 331 901 F: +420 296 331 935
Denmark (Freephone)
T: 80 30 80 02 F: 80 30 80 05
Eire (Freephone)
T: 1800 882 969 F: 1800 882 968
Finland
T: +358 20 762 8882 F: +358 20 762 8883
France
T: +33 1 39 30 49 00 F: +33 1 39 30 49 01
Germany (Freephone)
T: 08000 723300 F: 08000 723319
Greece
T: +30 210 9542 300 F: +30 210 9542 310
Hong Kong
T: +852 2960 7100 F: +852 2563 5276
Hungary
T: +36 1 456 30 40 F: +36 1 456 30 41
India – Bangalore
T: +91 80 2212 4928/9 F: +91 80 2207 5153
India – Chennai
T: +91 44 2220 1879 F: +91 44 2220 0716
India – Delhi
T: +91 11 4243 2222 F: +91 11 2610 0808

India – Mumbai
T: +91 22 4030 9165 F: +91 22 4040 9199
Italy
T: +39 010 60 60 1 F: +39 010 60 60 200
Korea
T: +82 2 3476 2121 F: +82 2 3476 1776
Mexico
T: +52 55 5687 63 90 F: +52 55 5687 62 28
T: +361 1 805 410 F: +361 1 805 445
Middle East/Gulf/North Africa
T: +31 30 602 55 44 F: +31 30 602 55 05
The Netherlands
T: 800 104 04 F: 800 101 90
Norway (Freephone)
T: 800 104 04 F: 800 101 90
Poland
T: +48 22 435 14 14 F: +48 22 435 14 10
Portugal
T: +351 1 381 25 40 F: +351 21 381 25 58
South Africa
T: +27 11 840 8600 F: +27 11 463 6077
South East Asia – Malaysia
T: +60 3 2283 3813 F: +60 3 2284 3813
South East Asia – Philippines
T: +63 2 687 6994 F: +63 2 687 3047
South East Asia – Singapore
T: +65 6418 8888 F: +65 6418 8899

South East Asia – Thailand
T: +66 2 6380 100 F: +66 2 6380 400
South East Europe
T: +30 210 9542 300 F: +30 210 9542 310
Spain
T: +34 901 11 12 15 F: +34 91 319 50 03
Sweden
T: +46 42 25 69 00 F: +46 42 25 69 69
Switzerland
T: +41 32 626 57 00 F: +41 32 626 57 01
Taiwan
T: +886 2 2747 7278 F: +886 2 2747 7270
Turkey
T: +90 216 464 36 66 F: +90 216 464 36 67
UK (Freephone)
T: 0844 800 4512 F: 0844 800 4513
Uruguay
T: +59 82 900 6212 F: +59 82 900 6212
Venezuela
T: +58 212 959 6275 F: +58 212 959 5328

Boston Scientific

Defining tomorrow, today.™

www.bostonscientific-international.com