

Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher,
NaviFlex Pancreatic Pusher

**Boston
Scientific**
Advancing science for life™

ORDERING INFO

STENT KITS

STRAIGHT STENT
— No Leading Barb

STRAIGHT STENT
— Leading Barb

PIGTAIL STENT
— No Leading Barb

PIGTAIL STENT
— Leading Barb

PUSHERS &
DELIVERY SYSTEM

Designed for Drainage, Flexibility, Durability
and Enhanced Visualization*

Part of a complete portfolio of devices for the treatment of pancreatico-biliary diseases.

*Radiopaque markers are not available on all stents.
Data claim does not apply to the 3Fr stent.

Learn more  The
National Pancreas
Foundation



ABOUT

FEATURES

DELIVERY

PACKAGING

GUIDEWIRES

Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher


ORDERING INFO

ABOUT

The Advanix Pancreatic Stent is indicated for pancreatic duct drainage and is designed to address **visualization** and **stent placement**.

In a systematic review involving 680 patients in 8 studies, pancreatitis was significantly reduced with the use of pancreatic duct stents from 19% in controls to 6%. The number needed for treatment in order to avoid a single episode of post-ERCP pancreatitis with pancreatic duct stent placement was 8.¹

¹ASGE Standards of Practice Committee, Complications of ERCP, GIE Journal, Volume 75, No. 3 : 2012



ABOUT

FEATURES

DELIVERY

PACKAGING

GUIDEWIRES

Boston
Scientific

Advanix™ Pancreatic Stent




and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher

ORDERING INFO

FEATURES: DRAINAGE

Drainage Holes:

- Facilitates flow into the drainage holes and flow into the duodenum
- Positioned every 5mm – at 90 degree intervals to support drainage and patency



* The Advanix 7Fr Pancreatic Stent inner diameter is 46% larger than the Cook 7Fr pancreatic stent.

*Testing conducted by Boston Scientific Research and Development.

Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher

ORDERING INFO

FEATURES: FLEXIBILITY AND DURABILITY

To aid in placement, the stent is constructed of a flexible, yet durable polymer.

“The nature of the material, the range of shapes and sizes and having radiopaque markers on the larger stent sizes will make this group of stents very useful for treating the pancreas.”

– **David L. Carr-Locke, M.D.**

Chief, Division of Digestive Diseases

Beth Israel Medical Center, New York, NY

< 2 of 3 >



ABOUT

FEATURES

DELIVERY

PACKAGING

GUIDEWIRES

Boston
Scientific


Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher

ORDERING INFO

FEATURES: ENHANCED VISUALIZATION*

- Endoscopic visible marker
- Radiopaque marker*
 - During procedure: Useful for proper stent placement
 - Post procedure: Designed to aid in visualization during follow-up fluoroscopic exam**



Fluoroscopic image from animal lab testing. Data on file, PDM 90895760.

*Radiopaque markers are not available on all stents. Data claim does not apply to the 3Fr stent.

**Bench test results may not necessarily be indicative of clinical performance.

Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher

ORDERING INFO

DELIVERY

The **NaviFlex Pancreatic Delivery System** and **NaviFlex Pancreatic Pushers** are designed to be pushable and durable — helping in stent placement, particularly in tortuous anatomy.



The Advanix Stent **NaviFlex RX Pancreatic Pusher** and **NaviFlex Pancreatic Pusher** are available separately or as a kit. Both offer short-wire (compatible with the RX Biliary System) and Long-Wire guidewire options for stent placement to accommodate physician preference.



ABOUT

FEATURES

DELIVERY

PACKAGING

GUIDEWIRES

Boston
Scientific


Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher

ORDERING INFO

PACKAGING


Boston Scientific is continually improving its products as well as packaging to address the clinical and procedural-efficiency needs of physicians, nurses and technicians.



Packing is color coded by French size, allowing nurses and technicians to quickly identify the desired size.

Example: Illustration of the 5Fr stent, straight with leading barb

Actual and working lengths — specifications help physicians with appropriate stent selection.



ABOUT

FEATURES

DELIVERY

PACKAGING

GUIDEWIRES

Boston Scientific

Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher

ORDERING INFO

GUIDEWIRES

Compatible guidewires include the Dreamwire™, Jagwire™ and Hydra Jagwire™ High Performance Guidewires.

“Having a four French pancreatic stent that can be placed over an .035" guidewire will be an advantage – one we haven't had until the Advanix Pancreatic Stent.”

– **David L. Carr-Locke, M.D.**

Chief, Division of Digestive Diseases

Beth Israel Medical Center, New York, NY

The Advanix 4Fr Pancreatic Stent is compatible with an .035in guidewire. Use of smaller caliber stents (3 or 4 French) have been shown to be associated with lower rates of duct injury than conventional polyethylene 5 French stents.²

²Freeman ML, Complications of Endoscopic Retrograde Cholangiopancreatography, Avoidance and Management. Gastrointest Endosc 2012; 05.001



ABOUT

FEATURES

DELIVERY

PACKAGING

GUIDEWIRES

Boston
Scientific

Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher

ORDERING INFO

ORDERING INFORMATION

STENT KITS

STRAIGHT STENT
— No Leading Barb

STRAIGHT STENT
— Leading Barb

PIGTAIL STENT
— No Leading Barb

PIGTAIL STENT
— Leading Barb

PUSHERS &
DELIVERY SYSTEM

| Code | Description | Diameter (Fr) | Length (cm) |
|-----------|--|---------------|-------------|
| M00537350 | Single Pigtail Stent – No Leading Barb | 4 | 3 |
| M00537360 | Single Pigtail Stent – No Leading Barb | 4 | 4 |
| M00537370 | Single Pigtail Stent – No Leading Barb | 4 | 5 |
| M00537380 | Single Pigtail Stent – No Leading Barb | 4 | 7 |
| M00537390 | Single Pigtail Stent – No Leading Barb | 5 | 3 |
| M00537400 | Single Pigtail Stent – No Leading Barb | 5 | 4 |
| M00537410 | Single Pigtail Stent – No Leading Barb | 5 | 5 |
| M00537420 | Single Pigtail Stent – No Leading Barb | 5 | 7 |
| M00537430 | Single Pigtail Stent – Leading Barb | 5 | 3 |
| M00537440 | Single Pigtail Stent – Leading Barb | 5 | 4 |
| M00537450 | Single Pigtail Stent – Leading Barb | 5 | 5 |
| M00537460 | Straight Stent – No Leading Barb | 4 | 3 |
| M00537470 | Straight Stent – No Leading Barb | 4 | 4 |
| M00537480 | Straight Stent – No Leading Barb | 4 | 5 |
| M00537490 | Straight Stent – No Leading Barb | 4 | 7 |
| M00537500 | Straight Stent – No Leading Barb | 5 | 3 |
| M00537510 | Straight Stent – No Leading Barb | 5 | 4 |
| M00537520 | Straight Stent – No Leading Barb | 5 | 5 |
| M00537530 | Straight Stent – No Leading Barb | 5 | 7 |
| M00537540 | Straight Stent – Leading Barb | 5 | 3 |
| M00537550 | Straight Stent – Leading Barb | 5 | 4 |
| M00537560 | Straight Stent – Leading Barb | 5 | 5 |

4 and 5Fr stents compatible with an .035" or smaller guidewire.

Scope compatibility: 7Fr stent requires a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.

5Fr RX Pusher with
RX Port and Guidewire



ABOUT

FEATURES

DELIVERY

PACKAGING

GUIDEWIRES

Boston
Scientific

Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher

ORDERING INFO

ORDERING INFORMATION

STENT KITS

STRAIGHT STENT
— No Leading Barb

STRAIGHT STENT
— Leading Barb

PIGTAIL STENT
— No Leading Barb

PIGTAIL STENT
— Leading Barb

PUSHERS &
DELIVERY SYSTEM

| Code | Diameter (Fr) | Length (cm) | Markers |
|-----------|---------------|-------------|----------------------|
| M00536000 | 3 | 2 | Trailing endo marker |
| M00536010 | 3 | 3 | |
| M00536020 | 3 | 4 | |
| M00536030 | 3 | 5 | |
| M00536040 | 3 | 6 | |
| M00536050 | 3 | 7 | |
| M00536060 | 3 | 8 | |
| M00536070 | 3 | 9 | |
| M00536080 | 3 | 10 | |
| M00536090 | 3 | 11 | |
| M00536100 | 3 | 12 | |
| M00536110 | 3 | 13 | |

| Code | Diameter (Fr) | Length (cm) | Markers |
|-----------|---------------|-------------|---|
| M00536120 | 4 | 2 | Trailing radiopaque marker |
| M00536130 | 4 | 3 | |
| M00536140 | 4 | 4 | Leading radiopaque marker with trailing endo marker |
| M00536150 | 4 | 5 | |
| M00536160 | 4 | 6 | |
| M00536170 | 4 | 7 | |
| M00536180 | 4 | 8 | |
| M00536190 | 4 | 9 | |
| M00536200 | 4 | 10 | |
| M00536210 | 4 | 11 | |
| M00536220 | 4 | 12 | |
| M00536230 | 4 | 13 | |

| Code | Diameter (Fr) | Length (cm) | Markers |
|-----------|---------------|-------------|---|
| M00536240 | 5 | 2 | Trailing radiopaque marker |
| M00536250 | 5 | 3 | |
| M00536260 | 5 | 4 | Leading radiopaque marker with trailing endo marker |
| M00536270 | 5 | 5 | |
| M00536280 | 5 | 6 | |
| M00536290 | 5 | 7 | |
| M00536300 | 5 | 8 | |
| M00536310 | 5 | 9 | |
| M00536320 | 5 | 10 | |
| M00536330 | 5 | 11 | |
| M00536340 | 5 | 12 | |
| M00536350 | 5 | 13 | |



3Fr stents compatible with an .021" or smaller guidewire. 4 and 5Fr stents compatible with an .035" or smaller guidewire. Scope compatibility: 7Fr stent requires a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



ABOUT

FEATURES

DELIVERY

PACKAGING

GUIDEWIRES

Boston Scientific

Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher

ORDERING INFO

ORDERING INFORMATION

STENT KITS

STRAIGHT STENT
— No Leading Barb

STRAIGHT STENT
— **Leading Barb**

PIGTAIL STENT
— No Leading Barb

PIGTAIL STENT
— Leading Barb

PUSHERS &
DELIVERY SYSTEM

| Code | Diameter (Fr) | Length (cm) | Markers |
|-----------|---------------|-------------|---|
| M00536360 | 5 | 2 | Trailing radiopaque marker |
| M00536370 | 5 | 3 | |
| M00536380 | 5 | 4 | Leading radiopaque marker with trailing endo marker |
| M00536390 | 5 | 5 | |
| M00536400 | 5 | 6 | |
| M00536410 | 5 | 7 | |
| M00536420 | 5 | 8 | |
| M00536430 | 5 | 9 | |
| M00536440 | 5 | 10 | |
| M00536450 | 5 | 11 | |
| M00536460 | 5 | 12 | |
| M00536470 | 5 | 13 | |

| Code | Diameter (Fr) | Length (cm) | Markers |
|-----------|---------------|-------------|---|
| M00536480 | 7 | 3 | Trailing endoscopic marker |
| M00536490 | 7 | 4 | Leading radiopaque marker with trailing endo marker |
| M00536500 | 7 | 5 | |
| M00536510 | 7 | 6 | |
| M00536520 | 7 | 7 | |
| M00536530 | 7 | 8 | |
| M00536540 | 7 | 9 | |
| M00536550 | 7 | 10 | |
| M00536560 | 7 | 11 | |
| M00536570 | 7 | 12 | |
| M00536580 | 7 | 13 | |
| M00536590 | 7 | 15 | |
| M00536600 | 7 | 18 | |

| Code | Diameter (Fr) | Length (cm) | Markers |
|-----------|---------------|-------------|---|
| M00536610 | 10 | 3 | Trailing endoscopic marker |
| M00536620 | 10 | 4 | Leading radiopaque marker with trailing endo marker |
| M00536630 | 10 | 5 | |
| M00536640 | 10 | 6 | |
| M00536650 | 10 | 7 | |
| M00536660 | 10 | 8 | |
| M00536670 | 10 | 9 | |
| M00536680 | 10 | 10 | |
| M00536690 | 10 | 11 | |
| M00536700 | 10 | 12 | |
| M00536710 | 10 | 13 | |
| M00536720 | 10 | 15 | |
| M00536730 | 10 | 18 | |



7 and 10Fr stents compatible with the NaviFlex Pancreatic Delivery System.
Scope compatibility: 7Fr stent requires a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



ABOUT

FEATURES

DELIVERY

PACKAGING

GUIDEWIRES

Boston
Scientific

Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher

ORDERING INFO

ORDERING INFORMATION

STENT KITS

STRAIGHT STENT
— No Leading Barb

STRAIGHT STENT
— Leading Barb

PIGTAIL STENT
— **No Leading Barb**

PIGTAIL STENT
— Leading Barb

PUSHERS &
DELIVERY SYSTEM

| Code | Diameter (Fr) | Length (cm) | Markers |
|-----------|---------------|-------------|----------------------|
| M00536740 | 3 | 2 | Trailing endo marker |
| M00536750 | 3 | 3 | |
| M00536760 | 3 | 4 | |
| M00536770 | 3 | 5 | |
| M00536780 | 3 | 6 | |
| M00536790 | 3 | 7 | |
| M00536800 | 3 | 8 | |
| M00536810 | 3 | 9 | |
| M00536820 | 3 | 10 | |
| M00536830 | 3 | 11 | |
| M00536840 | 3 | 12 | |
| M00536850 | 3 | 13 | |

| Code | Diameter (Fr) | Length (cm) | Markers |
|-----------|---------------|-------------|---|
| M00536860 | 4 | 2 | Trailing radiopaque marker |
| M00536870 | 4 | 3 | |
| M00536880 | 4 | 4 | Leading radiopaque marker with trailing endo marker |
| M00536890 | 4 | 5 | |
| M00536900 | 4 | 6 | |
| M00536910 | 4 | 7 | |
| M00536920 | 4 | 8 | |
| M00536930 | 4 | 9 | |
| M00536940 | 4 | 10 | |
| M00536950 | 4 | 11 | |
| M00536960 | 4 | 12 | |
| M00536970 | 4 | 13 | |

| Code | Diameter (Fr) | Length (cm) | Markers |
|-----------|---------------|-------------|---|
| M00536980 | 5 | 2 | Trailing radiopaque marker |
| M00536990 | 5 | 3 | |
| M00537000 | 5 | 4 | Leading radiopaque marker with trailing endo marker |
| M00537010 | 5 | 5 | |
| M00537020 | 5 | 6 | |
| M00537030 | 5 | 7 | |
| M00537040 | 5 | 8 | |
| M00537050 | 5 | 9 | |
| M00537060 | 5 | 10 | |
| M00537070 | 5 | 11 | |
| M00537080 | 5 | 12 | |
| M00537090 | 5 | 13 | |



3Fr stents compatible with an .021" or smaller guidewire. 4 and 5Fr stents compatible with an .035" or smaller guidewire. Scope compatibility: 7Fr stent requires a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



ABOUT

FEATURES

DELIVERY

PACKAGING

GUIDEWIRES

Boston Scientific

Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher

ORDERING INFO

ORDERING INFORMATION

STENT KITS

STRAIGHT STENT
— No Leading Barb

STRAIGHT STENT
— Leading Barb

PIGTAIL STENT
— No Leading Barb

PIGTAIL STENT
— Leading Barb

PUSHERS &
DELIVERY SYSTEM

| Code | Diameter (Fr) | Length (cm) | Markers |
|-----------|---------------|-------------|---|
| M00537100 | 5 | 2 | Trailing radiopaque marker |
| M00537110 | 5 | 3 | |
| M00537120 | 5 | 4 | Leading radiopaque marker with trailing endo marker |
| M00537130 | 5 | 5 | |
| M00537140 | 5 | 6 | |
| M00537150 | 5 | 7 | |
| M00537160 | 5 | 8 | |
| M00537170 | 5 | 9 | |
| M00537180 | 5 | 10 | |
| M00537190 | 5 | 11 | |
| M00537200 | 5 | 12 | |
| M00537210 | 5 | 13 | |

| Code | Diameter (Fr) | Length (cm) | Markers |
|-----------|---------------|-------------|---|
| M00537220 | 7 | 3 | Trailing endoscopic marker |
| M00537230 | 7 | 4 | Leading radiopaque marker with trailing endo marker |
| M00537240 | 7 | 5 | |
| M00537250 | 7 | 6 | |
| M00537260 | 7 | 7 | |
| M00537270 | 7 | 8 | |
| M00537280 | 7 | 9 | |
| M00537290 | 7 | 10 | |
| M00537300 | 7 | 11 | |
| M00537310 | 7 | 12 | |
| M00537320 | 7 | 13 | |
| M00537330 | 7 | 15 | |
| M00537340 | 7 | 18 | |



7 and 10Fr stents compatible with the NaviFlex Pancreatic Delivery System.
Scope compatibility: 7Fr stent requires a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



ABOUT

FEATURES

DELIVERY

PACKAGING

GUIDEWIRES

Boston Scientific

Advanix™ Pancreatic Stent

and NaviFlex™ RX Pancreatic Delivery System, NaviFlex RX Pancreatic Pusher, NaviFlex Pancreatic Pusher

ORDERING INFO

ORDERING INFORMATION

STENT KITS

STRAIGHT STENT
— No Leading Barb

STRAIGHT STENT
— Leading Barb

PIGTAIL STENT
— No Leading Barb

PIGTAIL STENT
— Leading Barb

PUSHERS & DELIVERY SYSTEM

| Code | Description | Diameter (Fr) | Length (cm) |
|-----------|--|---------------|-------------|
| M00535000 | Pancreatic Stent RX Pusher | 3 | n/a |
| M00535010 | Pancreatic Stent RX Pusher | 4 & 5 | n/a |
| M00535020 | Pancreatic Stent Pusher | 3 | n/a |
| M00535030 | Pancreatic Stent Pusher | 4 & 5 | n/a |
| M00535060 | NaviFlex RX Pancreatic Stent Delivery System | 7 | n/a |
| M00535080 | NaviFlex RX Pancreatic Stent Delivery System | 10 | n/a |



Delivery with Stent

4 and 5Fr stents compatible with an .035" or smaller guidewire.
Scope compatibility: 7Fr stent requires a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



ABOUT

FEATURES

DELIVERY

PACKAGING

GUIDEWIRES

Boston Scientific



Working with the National Pancreas Foundation

Boston Scientific is working with the National Pancreas Foundation to provide learning modules for patients. Using a mix of animations, physician testimonials, and downloadable information, the website is intended to provide easy-to-understand information for patients and enhance communication between a patient and their physician.

The website's learning modules could serve as an educational resource for health care professionals working with patients in a practice setting to help explain indications for the procedures, risks associated with the procedures, and the benefits of diagnostic and therapeutic pancreatic endoscopy, EUS FNA and ERCP.

www.pancreasfoundation.org

www.animatedpancreaspatient.com

Inclusion of information about the National Pancreas Foundation is not an endorsement by the National Pancreas Foundation of any Boston Scientific device.

**Boston
Scientific**

Advancing science for life™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com/endo-resources

Ordering Information
1.888.272.1001

©2014 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-246903-AB September 2014

Advanix, Dreamwire, Hydra Jagwire, Jagwire, and NaviFlex are unregistered or registered trademarks of Boston Scientific Corporation or its affiliates.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.