

SpyBite™ Max Biopsy Forceps

A design enhancement to the legacy SpyBite Biopsy Forceps have been shown to acquire **more than 2X tissue in an average bite**¹



**FRONT AND SIDE SERRATED
TEETH PROFILE**

**TWO ELONGATED FENESTRATION
HOLES**

INTERNAL SPIKE REMOVED



Performing biopsies under direct visualization using the SpyGlass System and SpyBite™ Biopsy Forceps (86% sensitivity) enables **faster, more accurate diagnosis of malignancies**²

1. Data on File
2. Shah et al. Performance of a fully disposable, digital, single-operator cholangiopancreatroscope. *Endoscopy* 2017; 49: 651–658.

SpyBite™ Max Biopsy Forceps

Redesigned Suite of Compatible Accessory Devices: Stricture Management

Getting More with SpyBite Max Biopsy Forceps

97%

TISSUE
ACQUISITION
RATE¹

Higher tissue acquisition rate can enable **more efficient** procedures

89%

OF
PHYSICIANS¹

Stated SpyBite Max acquires **more tissue** on an average bite than SpyBite

88%

OF
PHYSICIANS¹

Were **more confident** in their ability to acquire tissue from challenging scope positions

1. Data on File.

