


Rotatable Snare

Single-Use Polypectomy Snare

**Boston
Scientific**
Advancing science for life™

Rotatable snare is designed to provide **360 degree rotation** to help access difficult polyps.



Rotatable Snare

Single-Use Polypectomy Snare

FEATURES AND BENEFITS

Polyglide™ Sheath


- Predictable rotation both clockwise and counterclockwise
- Minimal friction when extending and retracting loop

Medium Stiff Braided Wire

- Designed to provide grip on polyps in difficult locations
- Control of wire

Hot / Cold Indication

- Allows clinical flexibility



Product Specifications

Loop Width - 13mm, 20mm

Sheath Diameter - 2.4mm

Working Length - 240cm

Working Channel - 2.8mm



Ordering Information

Order Number	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561821	Small Oval – Medium Stiff	13	2.4	2.8	240	Box 10
M00561831	Medium Oval – Medium Stiff	20	2.4	2.8	240	Box 10

All trademarks are the property of their respective owners.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

CAUTION: Federal (U.S.) law restricts this device to sale by or on the order of a physician.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Information for the use only in countries with applicable health authority product registrations.

**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com/gastro

©2016 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-393203-AA May 2016