

# AEON™

## RELOAD SELECTION

WIDEST  
RANGE  
OF  
ANVILS



STANDARD TIP



SHORT TIP



CURVED TIP

## RELOAD STAPLE HEIGHT

				
<b>GRAY</b>	<b>WHITE</b>	<b>ORANGE</b>	<b>PURPLE</b>	<b>BLACK</b>
OPEN 2.0	OPEN 2.5	OPEN 3.25	OPEN 4.0	OPEN 5.0
CLOSED 0.75	CLOSED 1.0	CLOSED 1.5	CLOSED 1.8	CLOSED 2.2
 2.0mm	 2.5mm	 3.25mm	 4.0mm	 5.0mm
 0.75mm	 1.0mm	 1.5mm	 1.8mm	 2.2mm





## 30mm RELOADS



## 45mm RELOADS



## 60mm RELOADS

Reload Type	Cartridge Lengths	Product Codes*	Open Height	Closed Height
Gray	30, 45, 60mm	AESR30G, AESR45G, AESR60G	2.0mm	0.75mm
Gray Curved Tip	30, 45mm	ASR30GC, AESC45G	2.0mm	0.75mm
Gray Short Tip	30mm	ASR30GS	2.0mm	0.75mm
White	30, 45, 60mm	AESR30W, AESR45W, AESR60W	2.5mm	1.0mm
White Curved Tip	30, 45, 60mm	ASR30WC, AESC45W, AESC60W	2.5mm	1.0mm
White Short Tip	30, 45, 60mm	ASR30WS, ASR45WS, ASR60WS	2.5mm	1.0mm
Orange	30, 45, 60mm	AESR30R, AESR45R, AESR60R	3.25mm	1.5mm
Orange Curved Tip	30, 45, 60mm	ASR30RC, AESC45R, AESC60R	3.25mm	1.5mm
Orange Short Tip	30, 45, 60mm	ASR30RS, ASR45RS, ASR60RS	3.25mm	1.5mm
Purple	30, 45, 60mm	AESR30P, AESR45P, AESR60P	4.0mm	1.8mm
Purple Curved Tip	30, 45, 60mm	ASR30PC, AESC45P, AESC60P	4.0mm	1.8mm
Purple Short Tip	30, 45, 60mm	ASR30PS, ASR45PS, ASR60PS	4.0mm	1.8mm
Black	45, 60mm	AESR45B, AESR60B	5.0mm	2.2mm
Black Curved Tip	45mm	ASR45BC	5.0mm	2.2mm



\*Some models may not be available in all geographies. © 2025 Lexington Medical, Inc. The third-party trademarks used herein are trademarks of their respective owners. MC-120723 Rev-D

Lexington Medical Inc.  
23 Crosby Drive  
Bedford, MA 01730 USA  
+1.617.209.9817



Adventa Ltd.  
Tower Business Center  
2nd Flr., Tower Street  
Swatar, BKR 4013 Malta