

Advancing science for life™

AMS 700[™]

Inflatable Penile Prosthesis



The AMS 700 Penile Prosthesis product line

The AMS 700 Prosthesis has an overall patient and partner satisfaction rate of **92%** and **96%**, respectively

Excellent girth expansion • Exceptional rigidity • Natural appearing flaccidity

Available with or without preconnected cylinders and pump for either an infrapubic or penoscrotal approach. Offering a combination of features focused on both surgeon and patient satisfaction.

DESIGNED FOR EASE OF IMPLANT/EASE OF USE

CONCEAL™ LOW PROFILE RESERVOIR

- Parylene coated to improve long-term durability of the device²
- Optimizes fit to offer a sleek, flexible solution

CYLINDER AND REAR TIP EXTENDER (RTE) DESIGN

- Reduced proximal profile designed to ease placement into proximal corpora even in difficult to dilate cases³
- Snap-fit RTEs designed to optimize fit and provide a secure, rigid connection
- Broad range of RTE lengths designed to custom fit patients (available in 0.5cm increments, ranging from 0.5–7.5cm)

MOMENTARY SQUEEZE (MS) PUMP

- One-touch button designed for deflation and ease of use
- Lock-out valve designed to resist auto inflation

INHIBIZONE™ ANTIBIOTIC SURFACE TREATMENT

- Only antibiotic-impregnated inflatable penile prosthesis on the market⁴
- Designed to create a zone of inhibition effective against the bacteria commonly associated with inflatable penile prosthesis infections⁵
- Provides a significant reduction in rate of revisions due to infection⁵

IMPORTANT SAFETY INFORMATION⁵

- Not for patients with urogenital infections or active skin infections in region of surgery
- InhibiZone Treatment version of this device is not for patients with known allergy or sensitivity to rifampin (rifampicin) or to minocycline, or other tetracyclines
- Will make latent natural or spontaneous erections, as well as other interventional treatment options impossible.
 Penile shortening, curvature or scarring may result
- · Patients with diabetes, spinal cord injuries or open sores may have an increased risk of infection associated with a prosthesis

AMS 700[™] Inflatable Penile Prosthesis

AMS 700

LGX INFLATABLE PENILE PROSTHESIS

AMS Men's Health exclusive: Cylinders designed to optimize **L**ength and **G**irth e**X**pansion with up to 20% elongation possible, depending on patients' anatomy¹

Product	InhibiZone [™] Item Number
LGX Preconnect MS 12cm, PS, IZ	72404250
LGX Preconnect MS 15cm, PS, IZ	72404251
LGX Preconnect MS 18cm, PS, IZ	72404252
LGX Preconnect MS 21cm, PS, IZ	72404253
LGX Preconnect MS 12cm, IP, IZ	72404255
LGX Preconnect MS 15cm, IP, IZ	72404256
LGX Preconnect MS 18cm, IP, IZ	72404257
LGX Preconnect MS 21cm, IP, IZ	72404258

AMS 700

CXR INFLATABLE PENILE PROSTHESIS

All the benefits of an IPP, including girth expansion, may be appropriate for patients who require smaller cylinders for those with scarred or sclerotic corpora¹

Product	InhibiZone Item Number
CXR Preconnect MS 10cm, PS, IZ	72404260
CXR Preconnect MS 12cm, PS, IZ	72404261
CXR Preconnect MS 14cm, PS, IZ	72404262
CXR Preconnect MS 16cm, PS, IZ	72404263
CXR Preconnect MS 18cm, PS, IZ	72404264
CXR Preconnect MS 10cm, IP, IZ	72404265
CXR Preconnect MS 12cm, IP, IZ	72404266
CXR Preconnect MS 14cm, IP, IZ	72404267
CXR Preconnect MS 16cm, IP, IZ	72404268
CXR Preconnect MS 18cm, IP, IZ	72404269

AMS 700

CX INFLATABLE PENILE PROSTHESIS

Controlled e**X**pansion cylinders that optimize girth and offer additional length options

Product	InhibiZone Item Number
CX Preconnect MS 12cm, PS, IZ	72404230
CX Preconnect MS 15cm, PS, IZ	72404231
CX Preconnect MS 18cm, PS, IZ	72404232
CX Preconnect MS 21cm, PS, IZ	72404233
CX Preconnect MS 24cm, PS, IZ	72404234*
CX Preconnect MS 12cm, IP, IZ	72404235
CX Preconnect MS 15cm, IP, IZ	72404236
CX Preconnect MS 18cm, IP, IZ	72404237
CX Preconnect MS 21cm, IP, IZ	72404238
CX Preconnect MS 24cm, IP, IZ	72404239*

^{*}Available as special order.

RESERVOIRS

Product	InhibiZone Item Number
Conceal™ Low Profile Reservoir	720185-01
Reservoir, 65mL, PC/IZ	72404155
Reservoir, 100mL, PC/IZ	72404156

^{*}All parts are also available without InhibiZone Antibiotic Surface Treatment.

All cited trademarks are the property of their respective owners.

CAUTION: Law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations. MH-401503-AA DEC 2016. Produced by Gosling.



www.bostonscientific.eu

© 2017 Boston Scientific Corporation or its affiliates. All rights reserved. DINURO2246EA

Montorsi F, Rigatti P, Carmignani G, et al. AMS three-piece inflatable implants for erectile dysfunction: a long-term multi-institutional study in 200 consecutive patients. Eur Urol. 2000 Jan;37(1):50-5.

^{3.} Data on file at Boston Scientific.

^{4.} FDA/PDP D970012-S065.

^{5.} AMS 700™ with MS Pump™ Penile Prosthesis Product Line Instructions for Use. American Medical Systems. 2013.