

Tactra[™]

Next-Generation Malleable Penile Prosthesis

The Natural Choice


A natural feeling malleable penile prosthesis constructed around a Nitinol core to optimise rigidity, durability and concealment.

- Proprietary dual-layer silicone construction for authentic, natural feel
- Dynamic Nitinol core provides optimal rigidity, durability and concealment
- Soft, rounded silicone distal tip designed to provide maximum patient comfort and satisfaction
- Simplified sizing to address the widest spectrum of patient anatomy
- Insertion-fit Rear Tip Extenders (RTEs) for a secure connection and stable foundation

• Ease of implantation through multiple surgical approaches

Has **4 times** the column strength required for intromission

The **only penile prosthesis** on the market with a Nitinol cable as its core



Simplified sizing

Three separate diameters trimmable to length intraoperatively

Diameter	Length	Product Number
9.5 mm	14 cm – 23 cm	720080-01
	Scientific Ex	\$1\$1\$1\$1\$1\$1\$1\$
11 mm	16 cm – 25 cm	720081-01
	Scientific Ex seconds	-11
13 mm	18 cm – 27 cm	720082-01
Scientific	2000000	-18-

Important safety information²

- Not for patients with urinary tract infections or active skin infections in region of surgery
- Will make latent natural or spontaneous erections, as well as other interventional treatment options, impossible
- Men with diabetes, spinal cord injuries or open sores may have increased risk of infection associated with a prosthesis

Please consult your sales representative for more information and ordering details.

- 1. Data on file with Boston Scientific
- 2. Tactra Malleable Penile Prosthesis Directions for Use. Boston Scientific. 2019
- 3. Al Ansari, A, Talib RA, Canguven O, et al. Axial penile rigidity influences patient and partner satisfaction after penile prosthesis implantation. *Arch Ital Urol Androl.* 2013 Sep 26;85(3):138-42.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for use only in countries with applicable health authority registrations. This material not intended for use in France. Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, and potential adverse events.

The Tactra™ Malleable Penile Prosthesis is intended for use in the treatment of erectile dysfunction (impotence) in adult males. Implanting a penile prosthesis will damage or destroy any remaining natural ability to have a spontaneous erection, as well as make other treatment options impossible.

Men with diabetes, spinal cord injuries, or skin infections may have an increased risk of infection. Implantation may result in penile shortening, curvature or scarring. Potential adverse events may include device malfunction/failure leading to additional surgery, device/tissue erosion, infection, and pain/soreness.

Product available in the European Economic Area (EEA) only. Please check availability with your local sales representative or customer service. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. Please check availability with your local sales representative or customer service.

All Images are the property of Boston Scientific. All trademarks are the property of their respective owners. The testing was performed by or on behalf of BSC. Data on file.

MH-586215-AA Produced by Gosling.



www.bostonscientific.eu

© 2019 Boston Scientific Corporation or its affiliates. All rights reserved. DINURO2407FA